

# TerraRack

Less Weight, Less Waste, 100% Recyclable

**With TerraRack™, you get a fresh clean rack of tips every time.**

TerraRack is a revolutionary new concept in pipette tip racks. They're half the weight of conventional racks, yet don't compromise on strength and are completely recyclable.

TerraRacks are single-use, disposable hinged-lid racks that won't flex or bend when loading single and multichannel pipettes – even when excessive force is applied. What's more, they compress easily for disposal and are largely composed of highly-recyclable PETE\*.

- Pre-sterilized (filter and non-filter)
- Individually wrapped
- 100% recyclable
- >50% less weight
- >50% less waste



\*PETE (polyethylene terephthalate) is perhaps the most commonly recycled plastic around the world. Widely used for water and soft-drink bottles, recycled PETE is in high demand by the beverage industry.

**METTLER TOLEDO**



## LTS TerraRack Tips for Rainin LTS® Pipettes

Order No.	Description	Max. Capacity	Price
<b>TerraRack LTS tips in PETE TerraRack with hinged lid</b>			
17014960	TR-L10S 960 LTS tips in 10 racks, presterilized	20 µL	\$ 66.25
17014961	TR-L10F 960 LTS tips in 10 racks, presterilized, filter	20 µL	102.50
17014962	TR-L250S 960 LTS tips in 10 racks, presterilized	200 µL	66.25
17014963	TR-L200F 960 LTS tips in 10 racks, presterilized, filter	200 µL	102.50
17014964	TR-L300S 768 LTS tips in 8 racks, presterilized	300 µL	66.25
17014965	TR-L300F 768 LTS tips in 8 racks, presterilized, filter	300 µL	102.50
17014966	TR-L1000S 768 LTS tips in 8 racks, presterilized	1000 µL	66.25
17014967	TR-L1000F 768 LTS tips in 8 racks, presterilized, filter	1000 µL	102.50
17014968	TR-L1200S 768 LTS tips in 8 racks, presterilized	1200 µL	66.25
17014969	TR-L1200F 768 LTS tips in 8 racks, presterilized, filter	1200 µL	102.50



## Universal TerraRack Tips

For Rainin®, Pipetman® & all brands with universal shaft design

Order No.	Description	Max. Capacity	Price
<b>TerraRack Tips in PETE TerraRack with hinged lid</b>			
17014970	TR-10S 960 tips in 10 racks, presterilized	10 µL	\$ 66.25
17014971	TR-10F 960 tips in 10 racks, presterilized, filter	10 µL	102.50
17014972	TR-10GS 960 ShaftGard tips in 10 racks, presterilized	10 µL	66.25
17014973	TR-10GF 960 ShaftGard tips in 10 racks, presterilized, filter	10 µL	102.50
17014974	TR-20F 960 tips in 10 racks, presterilized, filter	30 µL	102.50
17014975	TR-100F 960 tips in 10 racks, presterilized, filter	100 µL	102.50
17014976	TR-250S 960 tips in 10 racks, presterilized	250 µL	66.25
17014977	TR-200F 960 tips in 10 racks, presterilized, filter	200 µL	102.50
17014978	TR-300S 768 tips in 8 racks, presterilized	300 µL	66.25
17014979	TR-300F 768 tips in 8 racks, presterilized, filter	300 µL	102.50
17014980	TR-1000S 768 tips in 8 racks, presterilized	1000 µL	66.25
17014981	TR-1000F 768 tips in 8 racks, presterilized, filter	1000 µL	102.50

## Accessories

TerraBase™ for TerraRack Tips

Order No.	Description	Price
17014982	TR-SB TerraRack Base - short; 1 base fits 10 - 300 µL TR racks	\$ 15.00
17014983	TR-TB TerraRack Base - tall; 1 base fits 1000 µL TR racks	17.50

Rainin, Pipetting 360° and LTS are registered trademarks, and TerraRack and TerraBase are trademarks of Rainin Instrument, LLC.  
Pipetman is a registered trademark of Gilson, Inc.



### Rainin Instrument, LLC

7500 Edgewater Drive, Oakland, CA 94621  
Phone +1 510 564 1600

a METTLER TOLEDO Company



Subject to technical changes  
© 03/2015 Rainin Instrument, LLC  
Printed on demand from online content.